# The Revanesse® Family of Products

---

Revanesse® +  
Cross-Linked Hyaluronic Acid  
*with Lidocaine*

---

Instructions for Use  
Gebrauchsanweisung  
Mode d'emploi  
Istruzioni per l'uso  
Instrucciones de uso  
Instruções de Uso  
Инструкция по применению

Revitalize. Rejuvenate. ReNew®

## Table of Contents:

Symbols.....	4
Assembly of needle to syringe.....	6
Canada .....	8
English (en).....	8
Français (fr) .....	10
Brazil.....	12
Português (pt).....	12
Rest of the World.....	14
English (en).....	14
Deutsche (de).....	16
Français (fr) .....	18
Italiano (it) .....	20
Español (es) .....	22
Português (pt).....	24
Русский (ru) .....	26
EU.....	28
English (en).....	28
Deutsche (de).....	30
Français (fr) .....	32
Italiano (it) .....	34
Español (es).....	36
Português (pt).....	38

SYMBOL SYMBOL SYMBOLE SIMBOLO  
SIMBOLO SIMBOLOGIA СИМВОЛ:

Syringe fluid path sterilized using moist heat  
Sterilisation des Spritzenflüssigkeitspfades mittels feuchter Hitze  
Trajectoire du fluide dans la seringue stérilisée à la chaleur humide  
Percorso del fluido nella siringa sterilizzato a calore umido  
Canal del fluido de la jeringa esterilizado al vapor  
Esterilizado por vapor  
Путь протекающей по шприцу жидкости простерилизован с использованием влажного жара

Store between 2 and 25 °C  
Bei Temperaturen zwischen 2°C und 25 °C lagern  
Conserver entre 2 et 25 °C  
Conservare ad una temperatura compresa tra 2°C e 25°C  
Almacenar entre 2 y 25°C  
Armazenar entre 2 e 25 °C  
Храните при температуре от 2°C до 25°C

Read the Instructions before using the product  
Lesen Sie die Anweisungen, bevor Sie das Produkt verwenden  
Lire les instructions avant d'utiliser le produit  
Leggere le istruzioni prima di utilizzare il prodotto  
Lea las instrucciones antes de usar el producto  
Leia as instruções de uso antes de usar o produto  
Перед использованием изделия прочтите все инструкции

Do not use if the package has been damaged  
Verwenden Sie das Produkt nicht, wenn die Verpackung beschädigt ist  
Ne pas utiliser si l'emballage a été endommagé  
Non utilizzare se la confezione è danneggiata  
No usar si el paquete se encuentra dañado  
Não use se a embalagem estiver aberta ou danificada  
Не используйте, если упаковка повреждена

Do not reuse  
Das Produkt darf nicht wiederverwendet werden  
Ne pas réutiliser  
Evitare il riutilizzo  
No reutilizar  
Produto de uso único - proibido reprocessar - não reutilizar  
Не используйте повторно

Keep away from sunlight  
Vor Sonneneinstrahlung schützen  
Tenir à l'écart de la lumière du soleil  
Non esporre alla luce solare  
Mantener protegido de la luz solar  
Não expora luz solar  
Не допускайте попадания солнечного света

Sterilized by irradiation  
Strahlensterilisation  
Stérilisé par irradiation  
Sterilizzato tramite irradiazione  
Esterilizado por irradiación  
Esterilizado por irradiação  
Стерилизован облучением

Keep dry  
Trocken lagern  
Garder au sec  
Conservare in luogo secco  
Mantener seco  
Mantenha seco  
Храните в сухом месте

Lot number  
Chargennummer  
Numéro de lot  
Numero di lotto  
Número de LOTE  
Número do lote  
Номер партии

Expiration date  
Ablaufdatum  
Date d'expiration  
Data di scadenza  
Fecha de caducidad  
Data de validade  
Срок годности

Date of Manufacture  
Herstellungsdatum  
Date de fabrication  
Data di produzione  
Fecha de fabricación  
Data de Fabricação  
Дата изготовления

Manufacturer  
Hersteller  
Fabricant  
Produttore  
Fabricante  
Fabricante  
Производитель


SYRINGE  
SPRITZE  
SERINGUE  
SIRINGA  
JERINGA  
SERINGA  
ШПРИЦ

NEEDLE  
NADDEL  
AIGUILLE  
AGO  
AGUJA  
AGULHA  
ИГЛА

Emergo Europe B.V.  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

### ASSEMBLY OF NEEDLE TO SYRINGE:

1. Carefully unscrew the syringe end cap.
2. Firmly grip the wider part of the needle guard.
3. Press and rotate the needle guard by 90° (one quarter turn). This quarter turn is necessary to lock the needle into the syringe. Remove the needle guard.
4. Using the thumb and forefinger, tightly secure the glass syringe cylinder and the Luer-lock adapter.
5. Grip the needle guard with the opposite hand. To facilitate proper assembly, push and rotate firmly.

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
**COMPOSITION**

Cross-linked hyaluronic acid	25mg/ml
Lidocaine	0.3% w/w
In phosphate buffered saline	
[Cross linked with Butandiol-diglycidylether (BDDE)]	

**DESCRIPTION**

Revasense® + is a colorless, odorless, transparent and aqueous gel of synthetic origin. The gel is stored in a prefilled disposable syringe. Each box contains two 1ml syringes of Revasense® + along with two sterilized needles.

**APPLICATION RANGE / INDICATIONS**

Revasense® + is a cross linked hyaluronic acid gel that is indicated for the treatment of facial rhytides and lipotherapy by injection into the skin. Implantation life is dependent on depth and location of injection, and averages 6-12 months.

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- Injection-related reactions might occur, such as transient erythema, swelling, pain, itching, discoloration or tenderness at the injection site. These reactions may last for one week.
- Nodules or induration are also possible at the injection site.
- Poor product performance due to improper injection technique.
- Glabellar necrosis, abscess formation, granulomas and hypersensitivity have all been reported with injections of hyaluronic acid products. It is important for physicians to take these reactions into account on a case by case basis. Reactions thought to be of hypersensitivity in nature have been reported in less than one in every 1500 treatments. These have consisted of prolonged erythema, swelling and induration at the implant site.
- These reactions have started either shortly after injection or after a delay of 2 – 4 weeks and have been described as moderate, with an average duration of 2 weeks. Typically, this reaction is self-limiting and resolves spontaneously with time. However, it is imperative that patients with hypersensitivity type reactions contact their physician immediately for assessment. Patients with multiple allergic reactions should be excluded from the treatment.

**CONTRAINDICATIONS**

- Contains lidocaine and is contraindicated for patients with a history of allergies to such material.
- Do not inject Revasense® + into eye contours (into the eye circle or eyelids).
- Pregnant women, or women during lactation should not be treated with Revasense® +.
- Revasense® + is only intended for intradermal use and must not be injected into blood vessels. This may occlude and could cause an embolism.
- Patients who develop hypertrophic scarring should not be treated with Revasense® +.
- Contains trace amounts of gram positive bacterial proteins, and is contraindicated for patients with a history of allergies to such material.
- Never use Revasense® + in conjunction with a laser, intense pulsed light, chemical peeling or dermabrasion treatments.
- People under the age of 18 should not be treated with Revasense® +.
- Patients with acute and / or other inflammatory diseases of the skin should not be treated with Revasense® +.
- Patients with unrealistic expectations.
- Patients with auto-immune disorders or under immunotherapy.
- Patients with multiple severe allergies.
- Patients with acute or chronic skin disease in or near the injection sites.
- Coagulation defects or under anti-coagulation therapy.
- Patients with sensitivity to hyaluronic acid.
- It is imperative that patients with adverse inflammatory reactions that persist for more than one week report this immediately to their physician. These conditions should be treated as appropriate (ie: corticosteroids or antibiotics). All other types of adverse reactions should be reported directly to the authorized distributor of the Revasense® family of products and / or to ProLium Medical Technologies Inc. directly.

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**PRECAUTIONS**

- Revasense® + should not be injected into an area that already contains another filler product as there is no available clinical data on possible reactions.
- Revasense® + should not be injected into an area where there is a permanent filler or implant.
- Hyaluronic acid products have a known incompatibility with quaternary ammonium salts such as benzalkonium chloride. Please ensure that Revasense® + never comes into contact with this substance or medical instrumentation

that has come into contact with this substance.

- Revasense® + should never be used for breast enlargement, or for implantation into bone, tendon, ligament or muscle.
- Avoid touching the treated area for 12 hours after injection and avoid prolonged exposure to sunlight, UV, as well as extreme cold and heat.
- Until the initial swelling and redness have resolved, do not expose the treated area to intense heat (e.g. solarium and sunbathing) or extreme cold.
- If you have previously suffered from facial cold sores, there is a risk that the needle punctures could contribute to another eruption of cold sores.
- If you are using aspirin, non-steroidal anti-inflammatory medications, St. John's Wort or High doses of Vitamin E supplements prior to treatment or any similar medications be aware that these may increase bruising and bleeding at the injection site.
- Based on a toxicological risk assessment for lidocaine, patients should be limited to 20 mL per 60 kg (130 lbs) body mass per annum. The safety of injecting greater amounts has not been established.
- The safety for use in patients under 18 years or over 65 years has not been established.
- Patients who are visibly ill, with bacterial or viral infections, influenza, or active fever should not be treated until a resolution of their symptoms.
- In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.
- Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

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- Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction.
- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures.
- Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure.
- Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.
- The Revasense® family of products should not be used in areas that have high vascularity. Use in these areas such as the glabella and nose region has resulted

**SHELF LIFE & STORAGE**

Expiry is indicated on each individual package. Store between 2°-25° C, and protect from direct sun light and freezing.

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**MANUFACTURER**

ProLium Medical Technologies, Inc.  
138 Industrial Parkway North, Aurora, ON  
L4G 4C3, Canada

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**® FABRICANTE**

Prohlemium Medical Technologies, Inc.

138 Industrial Parkway North, Aurora, ON

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**COMPOSITION**

Cross-linked hyaluronic acid .....	25mg/ml
Lidocaine .....	0.3% w/w
In phosphate buffered saline	
[Cross linked with Butandiol-diglycidylether (BDDE)]	

**DESCRIPTION**

Revanesse® + is a colorless, odorless, transparent and aqueous gel of synthetic origin. The gel is stored in a pre-filled disposable syringe. Each box contains two 1ml syringes of Revanesse® + along with two sterilized needles.

**APPLICATION RANGE / INDICATIONS**

Revanesse® + is a cross linked hyaluronic acid gel that is indicated for the treatment of facial rhytides and lipatrophy by injection into the skin. Implantation life is dependent on depth and location of injection, and averages 6-12 months.

**ANTICIPATED SIDE EFFECTS**

Physicians must inform patients that with every injection of Revanesse® + there are potential adverse reactions that may be delayed or occur immediately after the injection. These include but are not limited to:

- Injection-related reactions might occur, such as transient erythema, swelling, pain, itching, discoloration or tenderness at the injection site. These reactions may last for one week.
- Nodules or induration are also possible at the injection site.
- Poor product performance due to improper injection technique.
- Gelabellar necrosis, abscess formation, granulomas and hypersensitivity have all been reported with injections of hyaluronic acid products. It is important for physicians to take these reactions into account on a case by case basis. Reactions thought to be of hypersensitivity in nature have been reported in less than one in every 1500 treatments. These have consisted of prolonged erythema, swelling and induration at the implant site.

These reactions have started either shortly after injection or after a delay of 2-4 weeks and have been described as mild to moderate, with an average duration of 2 weeks. Typically, this reaction is self-limiting and resolves spontaneously with time. However, it is imperative that patients with hypersensitivity type reactions contact their physician immediately for assessment. Patients with multiple allergic reactions should be excluded from the treatment.

**CONTRAINDICATIONS**

- Contains lidocaine and is contraindicated for patients with a history of allergies to such material.
  - Do not inject Revanesse® + into eye contours (into the eye circle or eyelids).
  - Pregnancy, or women during lactation should not be treated with Revanesse® +.
  - Revanesse® + is only intended for intradermal use and must not be injected into blood vessels. This may occlude and could cause an embolism.
  - Patients who develop hypertrophic scarring should not be treated with Revanesse® +.
  - Contains trace amounts of gram positive bacterial proteins, and is contraindicated for patients with a history of allergies to such material.
  - Never use Revanesse® + in conjunction with a laser, intense pulsed light, chemical peeling or dermabrasion treatments.
  - People under the age of 18 should not be treated with Revanesse® +.
  - Patients with acne and / or other inflammatory diseases of the skin should not be treated with Revanesse® +.
  - Patients with unattainable expectations.
  - Patients with auto-immune disorders or under immunotherapy.
  - Patients with multiple severe allergies.
  - Patients with acute or chronic skin disease in or near the injection sites.
  - Coagulation defects or under anti-coagulation therapy.
  - Patients with sensitivity to hyaluronic acid.
- It is imperative that patients with adverse inflammatory reactions that persist for more than one week report this immediately to their physician. These conditions should be treated as appropriate (i.e. corticosteroids or antibiotics). All other types of adverse reactions should be reported directly to the authorized distributor of the Revanesse® family of products and / or to Prolenium Medical Technologies Inc. directly.

**ADMINISTRATION & DOSAGE**

- Revanesse® + should only be injected by or under the direct supervision of qualified physicians who have been trained on the proper injection technique for filling facial wrinkles.
- Before patients are treated they should be informed of the indications of the device as well as its contraindications and potential undesirable side effects.
- The area to be treated must be thoroughly disinfected. Be sure to inject only under sterile conditions.
- Inject the product slowly and apply the least amount of pressure necessary.
- Revanesse® + and needles packaged with it are for single use only. Do not re-use. If re-used, there is a risk of infection or transmission of blood born diseases.
- Keep the product at room temperature for 30 minutes prior to injection.
- The amount of Revanesse® + that is injected will depend on the severity of the depression or the level of correction desired.
- If the skin turns a white color (blanching), the injection should be stopped immediately and the area should be massaged until the skin returns to its normal color.
- Before injecting, press on the plunger of the syringe until a small drop is visible at the tip of the needle.

**PRECAUTIONS**

- Revanesse® + should not be injected into an area that already contains another filler product as there is no available clinical data on possible reactions.
- Revanesse® + should not be injected into an area where there is a permanent filler or implant.
- Hyaluronic acid products have a known incompatibility with quaternary ammonium salts such as benzalkonium chloride. Please ensure that Revanesse® + never comes into contact with this substance or medical instrumentation that has come into contact with this substance.
- Revanesse® + should never be used for breast enlargement, or for implantation into bone, tendon, ligament

or muscle.

- Avoid touching the treated area for 12 hours after injection and avoid prolonged exposure to sunlight, UV, as well as extreme cold and heat.
- Until the initial swelling and redness have resolved, do not expose the treated area to intense heat (e.g. solarium and sunbathing) or extreme cold.
- If you have previously suffered from facial cold sores, there is a risk that the needle punctures could contribute to another eruption of cold sores.
- If you are using aspirin, non-steroidal anti-inflammatory medications, St. John's Wort or High doses of Vitamin E supplements prior to treatment or any similar medications be aware that these may increase bruising and bleeding at the injection site.
- Based on a toxicological risk assessment for lidocaine, patients should be limited to 20 mL per 60 kg (130 lbs) body mass per annum. The safety of injecting greater amounts has not been established.
- The safety for use in patients under 18 years or over 65 years has not been established.
- Patients who are visibly ill, with bacterial or viral infections, influenza, or active fever should not be treated until a resolution of their symptoms.
- In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy and at around the site of injection.
- Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

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- Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction.
- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures.
- Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure.
- Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.
- The Revanesse® family of products should not be used in areas that have high vascularity. Use in these areas such as the glabella and nose region has resulted in cases of vascular embolization and symptoms consistent with ocular vessel occlusion (i.e.: blindness).

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**NOTE:** The correct injection technique is crucial to treatment success & patient satisfaction. Revanesse® + should only be injected by a practitioner qualified according to local laws and standards.

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**MANUFACTURER**

Prolenium Medical Technologies, Inc.  
138 Industrial Parkway North, Aurora, ON  
L4G 4C3, Canada

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Prolium Medical Technologies, Inc.  
138 Industrial Parkway North, Aurora, ON  
L4G 4C3, Канада

**COMPOSITION**

Cross-linked hyaluronic acid	25mg/ml
Lidocaine	3mg/ml
In phosphate buffered saline [Cross linked with Butandiol-diglycidylether (BDDE)]	

**DESCRIPTION**

Revasense® is a colorless, odorless, transparent and aqueous gel of synthetic origin. The gel is stored in a pre-filled disposable syringe. Each box contains two 1ml syringes of Revasense® + along with two sterilized needles.

**APPLICATION RANGE / INDICATIONS**

Application: Superficial lines and wrinkles.

Medical Indications: The products are space-occupying tissue reconstructive materials composed of a hyaluronic acid gel that is indicated for restoration of volume lost from lipatrophy / lipodystrophy, and /or correction of contour deficiencies and anatomic deformities of either pathologic origin or after trauma, in soft tissue. Such intended patients are those desiring the correction of contour deficiencies and deformity in soft tissue, such as HIV-associated lipatrophy and lipodystrophy.

Cosmetic Indication: Revasense® + is indicated for the treatment of facial rhytides, volume restoration, lip augmentation, skin hydration and contouring of depressions by injection into tissue.

**ANTICIPATED SIDE EFFECTS**

Physicians must inform patients that with every injection of Revasense® + there are potential adverse reactions that may be delayed or occur immediately after the injection. These include but are not limited to:

- Injection-related reactions might occur, such as transient erythema, swelling, pain, itching, discoloration or tenderness at the injection site. These reactions may last for one week.
- Nodules or induration are also possible at the injection site.
- Poor product performance due to improper injection technique.
- Gelabular necrosis, abscess formation, granulomas and hypersensitivity have all been reported with injections of hyaluronic acid products. It is important for physicians to take these reactions into account on a case by case basis. Reactions to hyaluronic acid or to hypersensitivity in nature have been reported in less than one in every 1500 treated cases. These have consisted of prolonged erythema, swelling and induration at the implant site.
- These reactions have started either shortly after injection or after a delay of 2 – 4 weeks and have been described as mild or moderate, with an average duration of 2 weeks. Typically, this reaction is self-limiting and resolves spontaneously with time. However, it is imperative that patients with hypersensitivity type reactions contact their physician immediately for assessment. Patients with multiple allergic reactions should be excluded from the treatment.

**CONTRAINDICATIONS**

- Products are contraindicated in cases of hypersensitivity to all local amide-type anesthetics (not only to lidocaine as indicated in the package leaflet), as well as to all degrees of atrioventricular block, sinus node dysfunction, severe depression, myocardial conduction system, supraventricular tachycardias.
- Contains lidocaine and is contraindicated for patients with a history of allergies to such material.
- Do not inject Revasense® + into eye contours (into the eye circle or eyelids).
- Pregnant women, or women during lactation should not be treated with Revasense® +.
- Revasense® + is only intended for intradermal use and must not be injected into blood vessels. This may occlude and cause an embolism.
- Patients with contact dermatitis/hypersensitivity scarring should not be treated with Revasense® +.
- Contains trace amounts of gram positive bacterial proteins, and is contraindicated for patients with a history of allergies to such material.
- Never use Revasense® + in conjunction with a laser, intense pulsed light, chemical peeling or dermabrasion treatments.

- People under the age of 18 should not be treated with Revasense® +.
- Patients with acne and /or other inflammatory diseases of the skin should not be treated with Revasense® +.
- Patients with unattainable expectations.
- Patients with auto-immune disorders or under immunotherapy.
- Patients with multiple severe allergies.
- Patients with acute or chronic skin disease in or near the injection sites.
- Coagulation defects or under anti-coagulation therapy.
- Patients with sensitivity to hyaluronic acid.

It is imperative that patients with adverse inflammatory reactions that persist for more than one week report this immediately to their physician. These conditions should be treated as appropriate (i.e. corticosteroids or antibiotics). All other types of adverse reactions should be reported directly to the authorized distributor of the Revasense® family of products and /or to Prolium Medical Technologies Inc. directly.

**ADMINISTRATION & DOSAGE**

- Revasense® + should only be injected by or under the direct supervision of qualified physicians who have been trained on the proper injection technique for filling facial wrinkles.
- Before patients are treated they should be informed of the indications of the device as well as its contraindications and potential undesirable side effects.
- The area to be treated must be thoroughly disinfected. Be sure to inject only under sterile conditions.
- Inject the product slowly and apply the least amount of pressure necessary.
- Revasense® + and needles packaged with it are for single use only. Do not re-use. If re-used, there is a risk of infection or transmission of blood born diseases.
- Keep the product at room temperature for 30 minutes prior to injection.
- The amount of Revasense® + that is injected will depend on the severity of the depression or the level of correction desired.
- If the skin turns a white color (blanching), the injection should be stopped immediately and the area should be massaged until the skin returns to its normal color.
- Before injecting, press on the plunger of the syringe until a small drop is visible at the tip of the needle.

**PRECAUTIONS**

- If intravascular injections are made by mistake, the toxic effect will be noticed within 1-3 minutes (so the doctor and the patient should be alert in case of error).
- The effect of lidocaine may be reduced if infusion occurs in an inflammation or infection area.
- Special attention should be given to patients with partial or complete atrioventricular block because local anesthetics may suppress myocardial conduction in patients with advanced liver disease or severe renal

impairment, in patients with epilepsy, in patients with respiratory failure, in the elderly age, in patients with poor general health status, in patients receiving class III antiarrhythmics (eg amiodarone), who must be under close medical supervision including electrocardiogram due to the possible addition of cardiac effects, and finally in patients with acute porphyria.

- Lidocaine should be used with caution in patients receiving other topical anesthetics or agents with structural similarities to local amide-type anesthetics, e.g. certain antiarrhythmics, such as mexiletine and tocainide, since systemic toxic effects can be cumulative.

In addition to the direct anesthetic effect, local anesthetics may exert a very mild effect on the cognitive function and movement coordination and may temporarily affect somatometry and alertness. Depending on the dose, local anesthetics may have a very minor effect on mental function and may temporarily disturb movement and movement coordination.

- Revasense® + should not be injected into an area that already contains another filler product as there is no available clinical data on possible reactions.

- Revasense® + should not be injected into an area where there is a permanent filler or implant.
- Hyaluronic acid products have a known incompatibility with quaternary ammonium salts such as benzalkonium chloride. Please ensure that Revasense® + never comes into contact with this substance or medical instrumentation that has come into contact with this substance.

- Revasense® + should never be used for breast enlargement, or for implantation into bone, tendon, ligament or muscle.
- Avoid touching the treated area for 12 hours after injection and avoid prolonged exposure to sunlight, UV, as well as extreme cold and heat.

- Until the initial swelling and redness have resolved, do not expose the treated area to intense heat (e.g. solarium and sunbathing) or extreme cold.

- If you have previously suffered from facial cold sores, there is a risk that the needle punctures could contribute to another eruption of cold sores.

- If you are using aspirin, non-steroidal anti-inflammatory medications, St. John's Wort of High doses of Vitamin E supplements prior to treatment or any similar medications be aware that these may increase bruising and bleeding at the injection site.

- Based on a toxicological risk assessment for lidocaine, patients should be limited to a maximum of 20 mg per 60 kg (130 lbs) body mass per session. The safety of injecting greater amounts has not been established.

- This safety for use in patients under 18 years or over 65 years has not been established.
- Patients who are visibly ill, with bacterial or viral infections, influenza, or active fever should not be treated until a resolution of their symptoms.

- In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.

- Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

**WARNINGS**

Confirm that the seal on the box has not been broken and sterility has not been compromised. Confirm that the product has not expired. Product is for single use only; do not re-use. If re-used, there is a risk of infection or transmission of blood born diseases.

- Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction.

- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures.

- Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure.

- Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should receive intravascular injection occur.

- The Revasense® family of products should not be used in areas that have high vascularity. Use in these areas such as the glabella and nose region has resulted in cases of vascular embolization and symptoms consistent with ocular vessel occlusion (i.e.: blindness).

**SHELF LIFE & STORAGE**

Revasense® is indicated on each individual package. Store between 2°-25° C, and protect from direct sun light and freezing.

**NOTE:** The correct injection technique is crucial to treatment success & patient satisfaction. Revasense® + should only be injected by a practitioner qualified according to local laws and standards.

The graduation on the syringe is not precise and should be used as a guide only. The amount of material to be injected is best determined by visual and tactile assessment by the user.

**MANUFACTURER**

Prolium Medical Technologies, Inc.  
138 Industrial Parkway North, Aurora, ON  
L4G 4C3, Canada

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